

# Optafresh

Dextran-70 USP 0.2%  
Hypromellose 4000 BP 0.2%

## Sterile Eye Drops

**Composition:**

Each ml contains Dextran-70 USP 2 mg & Hypromellose 4000 BP 2 mg

**Description:**

It acts like natural tear, which provides soothing relief to dry, itchy, burning and irritated eyes. It helps eyes to stay moist, healthy and protects further irritation.

**Indications and usage:**

Optafresh used as artificial tear to prevent from damage to the cornea in patients with keratoconjunctivitis sicca and for ocular lubrication when Sjogren's syndrome affects lacrimal glands. It is also used for the temporary relief of burning and irritation due to dryness of the eye and for use as the protectant against further irritation, for temporary relief of discomfort due to expose to wind or sun.

**Dosage & Administration:**

Instill one drop into the affected eye(s) when required or as directed by the registered physician.

**Precautions:**

If allergic reactions are experienced, the use of drug should be discontinued and should consult with physician.!!

**Side-effects:**

No serious side effect is reported with the use of this drug, but in rare cases there may be allergic reactions.

**Storage :**

Store at 15°-20°C. Do not refrigerate or store at high temperatures. Keep container tightly closed. Keep out of reach of children.

**Packing :**

Each plastic dropper bottle contains 10 ml sterile Eye Drops.



Manufactured by  
**Apex Pharma Ltd.**  
(Eye Care Division)  
Shafipur, Gazipur