

Optafresh

Dextran 70 0.2% & Hypromellose 4000 0.2%

Sterile Eye Drops

Composition

Each ml contains Dextran 70 USP 2 mg and Hypromellose 4000 BP 2 mg.

Description

Optafresh acts like natural tear, which provides smoothing relief to dry, itchy, burning and irritated eyes. It helps eyes to stay moist, healthy and protects further irritation.

Indication

Optafresh used as artificial tear to prevent from damage to the cornea in patients with keratoconjunctivitis sicca and for ocular lubrication when Sjogren's syndrome affects lacrimal glands. It is also used for the temporary relief of burning and irritation due to dryness of the eye and for use as the protectant against further irritation, for temporary relief of discomfort due to exposure to wind or sun.

Dosage & Administration

Instill one drop into the affected eye when required or as directed by the registered physician.

Precaution

If allergic reactions are experienced, the use of drug should be discontinued and should consult with physician. Do not touch dropper tip to any surface to avoid contamination.

Side-effect

No serious side effect is reported with the use of this drug, but in rare cases there may be allergic reactions.

Storage

Store in cool, dry place and away from light. Do not refrigerate or store at high temperatures. Keep container tightly closed. Keep out of reach of children. Use within one month after the first opening.

Packaging

Each plastic dropper bottle contains 10 ml sterile Eye Drops.



A Syed Manzur Elahi Enterprise

Manufactured by
Apex Pharma Limited
Shafipur, Kaliakair, Gazipur, Bangladesh